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Software: Hystar and DataAnalysis (Bruker, Germany)

Ultimate 3000 RS Dionex-System (Dionex Softron, Germering, Germany) bestehend aus DGP-3600RS Pumpe, WPS-3000RS Rapid Separation Autosampler

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For the sum formula confirmation a flow injection analysis coupled to mass spectrometry (FIA-MS) was used. The LC-system consisted of an Ultimate 3000 RS Dionex system (Dionex Softron, Germering, Germany) with DGP-3600RS pump and WPS-3000RS Rapid Separation Autosampler. The high resolution mass spectrometer was a MicrOTOF-Q II (Bruker Daltonics, Bremen, Germany) with ESI interface (Bruker Daltonics, Bremen, Germany). The mobile phase of the LC-System consisted of acetonitrile : ammonium formate [10 mM] : formic acid XX:XX:0.1 (V/V/V) with a flow of 0.1 ml/min. The electrospray ion source was operated in positive ionization mode, scan range XXX – XXX m/z.

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ESI conditions: The capillary was set to ± 4.5 kV, the nebulizer was operated at 0.4 bar, the dry gas was set to 4 L/min at a temperature of 180 °C. Transfer voltages: The Funnel 1 and 2 were set to 200 Vpp. The hexapole RF voltage was set to 100 Vpp. The transfer time was adapted to 70 μ s. Mass calibration was done using a 20 mM sodium formate solution.

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MS/MS-experiments were accomplished in MRM-mode. The isolation width was set to 10 m/z and the collision energy to XX eV (siehe Ausdruck).